

Edge ONE™

A Powerful CO2 Laser for All Skin Conditions



Edge ONE™

Specifications

Laser Wavelength	10,600nm	
Laser Type	CO2 Laser	
Maximum Power	1 ~ 30W	
Aiming Beam	Diode Laser 655nm	
Surgical Handpiece	Laser Mode	CW, Single, Repeat, Group Pulse, ULTRA
	Power Output	1 ~ 30W, ULTRA
	Pulse duration	1 ~ 1000ms / ULTRA 0.08~1ms
	Pulse delay	1 ~ 1000ms
Fractional Handpiece	Frequency (ULTRA)	10~500Hz
	Power Output	1~30W, Edge Pulse
	Tip	120, 120(s), 350, 800
	Pulse Energy	1 ~ 300 mJ
	Density	~ 42.4 %
Dimensions (mm)	Scan Area	2 x 2 ~ 15 x 15 (mm)
	Scan Shapes	SQUARE, TRIANGLE, CIRCLE
Weight	997(H) x 240(W) x 282(L) 32kg	



Scan or visit to learn more.

Jeisys

Jeisys Medical Inc.

#401, Daeryung Techno Town 8th, 96, Gamasan-ro, Geumcheon-gu, Seoul, Republic of Korea
 Website : www.jeisys.com Tel : +82 2 2603 6417 Fax : +82 2 2603 6447 Email : sales@jeisys.com [L01A-C-B00]

Jeisys

Edge ONE™ / Surgical Mode

More Precise, More Safe

High-Peak Power

Ultra Pulse of Edge ONE specially makes the laser beam very accurate without damaging the surrounding tissue. It can also be applied to a very small lesions less than 1 millimeter.


INDICATION

Moles / Warts / Lentigo / Soft Fibroma
Syringoma / Actinic Cheilitis / Actinic Keratosis / More



Edge ONE™ / Fractional Mode


- Edge Pulse & Moving Mode
- Automatic Detection Tip
- Various Beam Pattern




INDICATION

Skin Tightening / Wrinkles, Fine lines / Pigmented Lesions & Lentigines
Scars (Acne Scar, Burn Scar, Hypertrophic Scar, Striae Distensae)
Fractional Ablative Skin Resurfacing / Large Pore / More

Various Beam Pattern



30W(Normal)




Edge ONE™ Edge Pulse



Edge Pulse delivers high peak power in a very short time and minimizes heat diffusion to reduce PIH.

Before treatment


3 months After 1 treatment



[Lentigo]

Before treatment


3 months After 1 treatment



[Viral Wart]

Before treatment

16 days After 1 treatment



[Burn Scar]

Before treatment

5 days After 1 treatment



[Rejuvenation]

Before treatment

4 weeks After 1 treatment



[Acne Scars]

Before treatment

4 weeks After the 2nd treatment



[Acne Scars]

Before treatment

4 weeks After the 2nd treatment



[Melasma]

Edge ONE™

A Powerful CO2 Laser for All Skin Conditions