

## LIPOCel™

### Specifications

Frequency	2MHz
Energy	Max. 150J/cm <sup>2</sup>
LCD	10.2 inch, touchscreen control panel display
Treatment Zone	30mm X 30mm
Focus Depth	8mm, 11mm, 13mm
Cooling	Chill cartridge for contact cooling(5°C)
Dimensions (mm)	1080(H) x 438(W) x 580(L)
Weight	80kg



Scan or visit to learn more.

# LIPOCel™

A New Paradigm for Body Shaping



**Jeisys** Jeisys Medical Inc.

#401, Daeryung Techno Town 8th, 96, Gamasan-ro, Geumcheon-gu, Seoul, Republic of Korea  
 Website : www.jeisy.com Tel : +82 2 2603 6417 Fax : +82 2 2603 6447 Email : sales@jeisy.com [H02A-C-800]



# LIPOcel™

" A New Paradigm for Body Shaping "

## [ 5°C ] Contact Cooling

Contact cooling protects epidermis and nerves from thermal damage, so it reduces pain dramatically and enables treatment with max energy up to **150J/cm²**.

### LIPOcel Cartridge



**Various depth cartridges** allow all area treatment

Cartridge Options  
 - No shot limit for 6 months  
 - No time limit for 3000 shots



5min.



15min.



10min.



8min.



6min.




6min.

### Contact Cooling

Safe Treatment - Temperature changes on epidermis



### Pain tolerance level changes with cooling




### Advantages

- High efficacy**  
150J(Max) energy induces both **Acute Effect** & **Chronic Apoptosis effect**
- Minimal pain**  
Contact cooling protects skin from thermal damage so reduces pain
- No downtime**  
Undamaged skin allows normal routine life
- Fast treatment**  
High energy & Contact cooling pain control enable 15min fast treatment
- Minimal running cost**  
Unlimited shot cartridge  
No consumable cost for 6 months

### Most Profitable Device

Smart ROI - The most profitable device for your successful business



**LIPOcel**  
With Unlimited Cartridge, no consumable cost occurs in 6 months.

**Competitors:**  
Consumable cost occurs for every single shot and treatment.

**PROFITABILITY** of LIPOcel is a lot higher than competitors



" Current study showed increasing number of patient seeking non-invasive procedure for body sculpting indication. After trying and comparing so many other energy-based treatment I noticed that **LIPOcel™** using HIFU technology is so far delivering the **fast, painless and most importantly effective result within short period of time**. It has scientific evidence on acute effect to tighten the skin and chronic effect to lysis subcutaneous fat cell to achieve the result my patient desire "

- Dr. Ungku Shahrin, MD, Medical Aesthetician, Malaysia -

### Lunch time treatment for busy life style

A New paradigm of body shaping treatment inspired by doctors

**15 min.**  
full abdomen treatment

- ✓ Easy & Comfort
- ✓ Triple Energy
- ✓ Lunch time
- ✓ Good result
- ✓ No Pain

